



Lagdu Singh Charitable Trust's (Regd.)

## THAKUR COLLEGE OF ENGINEERING & TECHNOLOGY

(Approved by AICTE, Govt. of Maharashtra & Affiliated to University of Mumbai\*)  
Institute Accredited by National Assessment and Accreditation Council (NAAC), Bangalore#

Programmes Accredited by National Board of Accreditation (NBA), New Delhi\*\*  
Among Top 200 Colleges in the Country where Ranked 193<sup>rd</sup> in NIRF India Ranking 2019  
in Engineering College category

\*Permanent Affiliated UG Programmes : • Computer Engineering • Electronics & Telecommunication Engineering • Information Technology (w.e.f. A.Y. 2015-16)  
• Electronics Engineering (w.e.f. A.Y. 2017-18)  
\*\*3rd cycle NBA Accredited UG Programmes : • Computer Engineering • Electronics & Telecommunication Engineering • Information Technology (3 years w.e.f. 01-07-2019)  
\*\*2nd cycle NBA Accredited UG Programmes : • Electronics Engineering (3 years w.e.f. 01-07-2019)  
# 1st cycle of NAAC Accreditation : • "A" Grade for 5 years (w.e.f. 30-10-2017)

A - Block, Thakur Educational Campus,  
Shyamnarayan Thakur Marg, Thakur Village,  
Kandivali (East), Mumbai - 400 101.

Tel.: 6730 8000 / 8106 / 8107

Fax : 2846 1890

Email : [tcet@thakureducation.org](mailto:tcet@thakureducation.org)

Website : [www.tcetmumbai.in](http://www.tcetmumbai.in) • [www.thakureducation.org](http://www.thakureducation.org)



ISO 9001:2015 Certified

Subject :- Chemistry

Experiment / Tutorial / Assignment No. :- 6

Page :- 1

Date :- / /

# DETERMINATION OF MOLAR EXTINCTION COEFFICIENT





- Aim: To determine molar extinction coefficient for known and unknown solution of potassium permanganate solution at  $\lambda_{max}$  using Colorimeter.

- Objectives:

After performing the practical, the learner will be able to

PRO 1: Understand the theory behind coloured  $KMnO_4$  solution.

PRO 2: Understand the uses of Colorimeter and Glassware.

PRO 3: Determine the absorbance of solution.

PRO 4: Determine Molar Extinction coefficient.

- Apparatus:

→ Colorimeter

→ 250 mL & 100 mL volumetric flask.

→ 0.01 N  $KMnO_4$  solution.

## Determination of Molar Extinction Coefficient at $\lambda_{max}$ using Colorimeter.

### Observation:

#### PRO 1:

a] The relation between absorbance and transmittance is given by

$$A = -\log_{10} T \text{ equation}$$

b] It has been observed that with increasing concentration, absorbance of solution is increasing.

#### PRO 2:

a] Light source used in colorimeter is tungsten filament lamp having polychromatic light whereas most commonly used detector is Photovoltaic cell.

b] Cuvette or sample holder used in colorimetric analysis can be made up glass because it does not affect or absorb in visible region.

#### PRO 3:

#### Observation table:

Obs. No.	Concentration of $\text{KMnO}_4$ sol <sup>n</sup> in $\text{m/dm}^3$	Absorbance of Soln. "A"	Molar Extinction Coefficient $\epsilon = A/bc$
1	$4 \times 10^{-5}$	0.17	$0.0425 \times 10^5$
2	$8 \times 10^{-5}$	0.38	$0.0475 \times 10^5$
3	$12 \times 10^{-5}$	0.57	$0.0475 \times 10^5$
4	$16 \times 10^{-5}$	0.68	$0.0418 \times 10^5$
5	Unknown sol <sup>n</sup> .	0.55	$0.0475 \times 10^5$





#### PRO 4:

Molar Extinction Co-efficient:

$$\epsilon = \frac{A}{bc}$$

$$\text{Unit} = \text{Lit/mol cm}$$

(1)  $\epsilon$  for concentration (c)  $4 \times 10^{-5} \text{ mol/dm}^3$  and Absorbance (A) 0.17

$$\epsilon = \frac{A}{bc} = \frac{0.17}{4 \times 10^{-5}} = 0.0425 \times 10^5 \text{ Lit/mol cm}$$

(2)  $\epsilon$  for concentration (c)  $8 \times 10^{-5} \text{ mol/dm}^3$  and Absorbance (A) 0.38

$$\epsilon = \frac{A}{bc} = \frac{0.38}{8 \times 10^{-5}} = 0.0475 \times 10^5 \text{ Lit/mol cm}$$

(3)  $\epsilon$  for concentration (c)  $12 \times 10^{-5} \text{ mol/dm}^3$  and Absorbance (A) 0.57

$$\epsilon = \frac{A}{bc} = \frac{0.57}{12 \times 10^{-5}} = 0.0475 \times 10^5 \text{ Lit/mol cm}$$

(4)  $\epsilon$  for concentration (c)  $16 \times 10^{-5} \text{ mol/dm}^3$  and Absorbance (A) 0.68

$$\epsilon = \frac{A}{bc} = \frac{0.68}{16 \times 10^{-5}} = 0.0418 \times 10^5 \text{ Lit/mol cm}$$

(5)  $\epsilon$ ,  $0.0475 \times 10^5$  for concentration Unknown (c), by Graph and Absorbance (A) 0.55

$$\epsilon = \frac{A}{bc} \Rightarrow 0.0475 \times 10^5 = \frac{0.55}{c}$$

$$\Rightarrow c = \frac{0.55}{0.0475 \times 10^5}$$

$$\Rightarrow c = 11.57 \times 10^{-5} \text{ mol/dm}^3$$



- Result & Discussions:

PRO 1:  $\text{KMnO}_4$  Solution when Incident with Ray of White Light, Reflects Purple Color while Absorbing Everything Else, Hence to Human Eye the Solution Looks Purple. This color is due to Transitions of charge within the Compound. As on Dilution the Color changes from Purple Tint to Pink.

PRO 2: Colorimeter is an Analytical Instrument based on Beer-Lambert's Law & is used to find the Concentration of a Solution by Measurement of its Relative Absorption of Light. It consists of a Straight Arrangement of Point Source of visible Light passing through Alterable Filters producing beams of Light of Different wavelength which are incident on a cleansed Cuvette containing the solution and a Detector setup to Analyze the Absorption.

PRO 3: The Absorption of Light by  $\text{KMnO}_4$  Solution is obtained by Giving Colorimeter a Blank, i.e. Distilled water & then the solution. This procedure is repeated to find Absorption of all Possible Wavelengths in Colorimeter. Precise Handling of the Cuvette is Required as Fingerprints while Holding can Alter the Readings. The wavelength at which the Absorbance is found to be Highest, is called  $\lambda_{\text{max}}$  i.e. Peak Wavelength.

PRO 4: Absorbance of Solutions with different Concentrations along with the Unknown is Calculated using,

Colorimeter at Molar Extinction. Peak Wavelength.  
By obtaining  $\lambda_{max}$ . First, the Procedure is simplified considerably. Molar Extinction Coefficient -  $\epsilon$  is calculated for solution of Each Concentration & Compared & Recurring value is considered ideal. The Concentration of Unknown is then Calculated using Graph.

### • Conclusion :

The Molar Extinction Coefficient -  $\epsilon$  for  $KMnO_4$  solution is  $0.0475 \times 10^5 \text{ Lit/mol/cm}$ . The Concentration of Unknown Solution is  $11.6 \times 10^{-5} \text{ m/dm}^3$ .

### • Precautions:

- While performing the practical, the learner must:
1. Ensure that chemicals must be handled carefully.
  2. Ensure selection of proper filter.
  3. Ensure use of blank filter or Wavelength changed.
  4. Safe handling of Cuvette.

### • Quiz:

1] What is equation of Beer's Lambert Law? Explain the term involved in the Equation:

⇒ The Beer-Lambert's Law states that,

Absorbance of Light of Color Compounds is Directly proportional to Concentration & Path Length.

$$A \propto b.C$$



$$A = \epsilon \cdot b \cdot C$$

where,  $A$  = Absorbance of Light for the given Compound  
 $b$  = Path Length covered by the Light  
 $C$  = Concentration of given Solution.

2] How will you adjust the Wavelength?

⇒ The Wavelength of Light can be Adjusted using a Filter by Passing a Beam of Light through it. A Filter only allows a certain wavelength to pass through it and blocks the other wavelength. A standard Colorimeter has a Setup of 8 Alterable Filters to obtain visible Light Rays of Different Wavelength.

3] Path length of Cuvette used in this experiment is:?

⇒ Path length of Cuvette used in this experiment is 10mm.

ISO 9001:2015 Certified

NBA and NAAC Accredited

4] Is this technique useful for measuring concentration of colourless solution?

⇒ As Colorimeter is used to determine concentration of coloured solution it is not possible to measure concentration of colourless solution using this technique.

Objective	PRO	PRO	PRO	PRO	
	1	2	3	4	Total
Weight					Score
Points					
Score					
Earned Points (EP) =	Marks in 100 = EP * 20				
Total score / 20 =	=				



Subject : \_\_\_\_\_

Topic : Beer-Lambert Law

Y-axis:  $\text{cm} = 0.1$

