



Laxmi Singh Charitable Trust's (Regd.)

THAKUR COLLEGE OF ENGINEERING & TECHNOLOGY

(Approved by AICTE, Govt. of Maharashtra & Affiliated to University of Mumbai*)
Institute Accredited by National Assessment and Accreditation Council (NAAC), Bangalore#
Programmes Accredited by National Board of Accreditation (NBA), New Delhi**
Among Top 200 Colleges in the Country where Ranked 193rd in NIRF India Ranking 2019
in Engineering College category

*Permanent Affiliated UG Programmes : • Computer Engineering • Electronics & Telecommunication Engineering • Information Technology (w.e.f. A.Y. 2015-16)
• Electronics Engineering (w.e.f. A.Y. 2017-18)
**3rd cycle NBA Accredited UG Programmes : • Computer Engineering • Electronics & Telecommunication Engineering • Information Technology (3 years w.e.f. 01-07-2019)
**2nd cycle NBA Accredited UG Programmes : • Electronics Engineering (3 years w.e.f. 01-07-2019)
1st cycle of NAAC Accreditation : • "A" Grade for 5 years (w.e.f. 30-10-2017)

A - Block, Thakur Educational Campus,
Shyamnarayan Thakur Marg, Thakur Village,
Kandivali (East), Mumbai - 400 101.

Tel.: 6730 8000 / 8106 / 8107

Fax : 2846 1890

Email : tcet@thakureducation.org

Website : www.tcetmumbai.in • www.thakureducation.org



ISO 9001:2015 Certified

Subject :- Chemistry

Experiment / Tutorial / Assignment No. :- 7

Page :- 1

Date :- 13/1/2023

SAPONIFICATION VALUE OF OIL



- Aim:

To determine the Saponification value of Oil using simulator.

- Objectives:

After performing the practical, the learner will be able to:

PRO1: Understand the concept of Saponification Value.

PRO2: Use the Animation provided in simulator for carrying out experiments on simulator.



PRO3: Use the Simulator for carrying out experiment for Blank and Test sample.

PRO4: Calculate Saponification value of all Test Samples.

- Reactions:

ISO 9001:2015 Certified

NBA and NAAC Accredited



- Apparatus:

- Burette
- Pipette
- Conical flask
- Measuring Cylinder

- Reagents:

- KOH 0.5N KOH
- 0.1N HCl 0.5N HCl
- Phenolphthalein.
- Fat Solvent [ethanol/ether mixture]

PART-I

(For Blank solution)

• Observations:

- Solution in Burette: 0.5 N HCL
- Solution in Conical Flask: Fat solvent + KOH
- Indicator: Phenolphthalein
- End point: Pink to Colourless
- Pilot Reading: 19.0 ml to 20.0 ml

• Observation Table:

Obs. No.	Initial Reading(ml)	Final Reading(ml)	Difference(ml)
1	0.0	20.0	20.0
2	0.0	20.0	20.0
3	0.0	20.0	20.0


Constant Burette Reading = 20.0 ml

Estimation of Saponification Value.

(For Oil solution)

• Observations:

- Solution in Burette: 0.5 N HCL
- Solution in Conical Flask: Test Oil Sample + Fat solvent + KOH
- Indicator: Phenolphthalein
- End point: Pink to Colourless
- Pilot Reading: 10.0 ml to ~~11.0~~ ml
11.0



ISO 9001:2015 Certified
NBA and NAAC Accredited

• Observation Table:

• Test Sample 1: Coconut Oil

Obs. No.	Initial Reading (ml)	Final Reading (ml)	Difference (ml)
1	0.0	11.0	11.0
2	0.0	11.0	11.0
3	0.0	11.0	11.0

• Constant Burette Reading = 11.0 ml

• Calculations:

$$\text{Sap value} = \frac{(\text{Blank reading} - \text{Back reading}) \times \text{Normality of KOH} \times 56}{\text{Weight of oil taken}}$$

$$= \frac{(20 - 11) \times 0.5 \times 56}{1000}$$

$$= 0.252 \text{ g} = \underline{252 \text{ mg}}$$

• Test Sample 2: Sunflower Oil

Obs. No.	Initial Reading (ml)	Final Reading (ml)	Difference (ml)
1	0.0	13.0	13.0
2	0.0	13.0	13.0
3	0.0	13.0	13.0

• Constant Burette Reading = 13.0 ml

• Calculations:

$$\text{Sap value} = \frac{(\text{Blank reading} - \text{Back reading}) \times \text{Normality of KOH} \times 56}{\text{Weight of oil taken}}$$


$$= \frac{(20 - 13) \times 0.5 \times 56}{1000}$$

$$= 0.196 \text{ g} = \underline{196 \text{ mg}}$$

Subject :- Chemistry

Experiment / Tutorial / Assignment No. :- 7

Page :- 4



ISO 9001:2015 Certified
NBA and NAAC Accredited

• Test Sample 3: Butter

Obs. No.	Initial Reading (ml)	Final Reading (ml)	Difference (ml)
1	0.0	12.5	12.5
2	0.0	12.5	12.5
3	0.0	12.5	12.5

• Constant Burette Reading = 12.5 ml

• Calculations:

$$\text{Sap value} = \frac{(\text{Blank reading} - \text{Back reading}) \times \text{Normality of KOH} \times 56}{\text{Weight of oil taken}}$$

$$= \frac{(20 - 12.5) \times 0.5 \times 56}{1000}$$

$$= 0.210 \text{ g} = \underline{210 \text{ mg}}$$

• Test Sample 4: Rice Bran Oil

Obs. No.	Initial Reading (ml)	Final Reading (ml)	Difference (ml)
1	0.0	13.3	13.3
2	0.0	13.3	13.3
3	0.0	13.3	13.3

• Constant Burette Reading = 13.3 ml

• Calculations:

$$\text{Sap value} = \frac{(\text{Blank reading} - \text{Back reading}) \times \text{Normality of KOH} \times 56}{\text{Weight of oil taken}}$$


$$= \frac{(20 - 13.3) \times 0.5 \times 56}{1000}$$

$$= 0.1876 \text{ g} = \underline{187.6 \text{ mg}}$$

Subject :- Chemistry

Experiment / Tutorial / Assignment No. :- 7

Page :- 5



ISO 9001:2015 Certified
NBA and NAAC Accredited

• Test Sample 5 : Castor oil

Obs. No.	Initial Reading (ml)	Final Reading (ml)	Difference (ml)
1	0.0	13.5	13.5
2	0.0	13.5	13.5
3	0.0	13.5	13.5

• Constant Burette Reading = 13.5 ml

• Calculations:

$$\text{Sap value} = \frac{(\text{Blank reading} - \text{Back reading}) \times \text{Normality of KOH} \times 56}{\text{Weight of oil taken}}$$

$$= \frac{(20 - 13.5) \times 0.5 \times 56}{1000}$$

$$= 0.182 \text{ g} = \underline{182 \text{ mg}}$$

• Result & Discussions:

PRO 1: Fats (triglycerides) upon Alkaline Hydrolysis (either with KOH or NaOH) yield glycerol and Potassium or Sodium salts of fatty acids (soap). This process is called as Saponification. The Saponification Number is the number of milligrams of Potassium Hydroxide - KOH required to neutralize the fatty acids resulting from the complete hydrolysis of 1 g fat. As KOH is used in Excess some of it remains Unreacted. This Amount can be found using Redox Titration to obtain Saponification Number.

PRO 2: The Procedure of Experiment is explained through Animation. Adding Ether/Ester Mixture Solution i.e. Fat Solvent in the oil sample bit by bit Results in Complete Dissolution of the fats present in the Oil. Reflux Condensers are used for Cooling the Sample Vapors to Turn them back into Liquid state, thus Reducing the Loss of Sample.

PRO 3: The Saponification Experiment is Vast & requires a lot of Time to perform for its Tedious Procedure, Hence it is not advisable to be conducted in a Single Setting. Hence we perform this Experiment using a Simulator. Use of simulator Not only saves time but also Nullifies every possible Human Error that can happen while performing the Experiments Manually, Resulting in Accurate Readings.

PRO 4: The Saponification Value is calculated as the Number of mg of KOH consumed by 1g of Fat. The Blank Solution is made using just the Fat Solvent & KOH & since No Hydrolysis Take place, KOH solution remains Unreacted, Causing its Complete Titration Against HCl. Hence the Saponification value of Blank is Always Greater than Any oil sample.

Conclusion :

The Saponification Value of Any Triglyceride (Fat) can be found using its Ability to Hydrolyse with Alkali. The higher the Saponification number, the more capable the Oil is in making Soap.

Precautions:

While performing the practical, the learner must:

1. Ensure that glassware is Clean.
2. Ensure Burette is holding properly to burette stand.
3. Ensure Safe handling of all glassware used.
4. Ensure safe handling of Alcoholic KOH and Oil solution on water bath.

Quiz:

- 1] What is the role of Phenolphthalein Indicator?
⇒ Phenolphthalein is a Visual Indicator which gives Pink Color to an Alkaline Solution & stays Colourless for Acidic Solutions. As this change of colour can be Observed while

changing the pH, Phenolphthalein is Used as a Indicator in Titrations.

2] Why heating of oil solution is required?

⇒ Heating of Solution has a Big Role in Saponification Reaction as when the Reaction Takes place at Temperature Higher than 120° , it helps in Easier & Better Saponification Also Increasing its Stability.

3] Why water bath is needed for heating?

⇒ In the Saponification Experiment the Fat Solutions we use are Highly Flammable Chemicals. Hence Water baths are used as Preferred source of Heating these solutions to prevent Direct Ignition.

4] Name the type of titration involved in Saponification of Oil.

⇒ Acid-Base Titration using Phenolphthalein Indicator is used while Performing Saponification of Oil.

Objective	PRO	PRO	PRO	PRO	
	1	2	3	4	Total
Weight	20	20	20	20	Score
Points					
Score					
Earned Points (EP) =		Marks in 100 =			
Total Score / 80 =		EP * 20			